

CABLE TANGO RS232 INTERFACE FOR CASE 8000 BOX 16...

SUNTECH

Ref. [DSU230748](#)

Cable TANGO RS232 interface for case 8000 box 16...

Description

Cable TANGO RS232 interface for case 8000 box 16...

Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	3700446012410	HS Code	9018110000
Warranty (months)	3	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	20	Packed unit / Height (cm)	5
Packed unit / Length (cm)	22	Packed unit / Weight	0.12

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---