

AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR SPENGLER CORAL

Ref. SPG345



CE
0459



08/10/2025 00:51:54 – For information purposes only, can be subject to change

AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR SPENGLER CORAL



Ref. **SPG345**

SPENGLER precision in electronic version

Description

MEASUREMENT
ACCURACY• SPENGLER MEASUREMENT•
WHO risk indicator

EASY &
CONVENIENT• Fully automatic• 2

buttons• LCD display
• 120 time-stamped memories•
Fast and quiet

ADDITIONAL INFORMATION•
Wrist circumference 13.5-21 cm
• Comes with 2 LR03 AAA batteries and a storage case

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446041595	HS Code	90189010
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	10	Packed unit / Height (cm)	10
Packed unit / Length (cm)	9	Packed unit / Weight	0.27
Bare Unit / Weight (kg)	0.228		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-25	Max temperature	70
Minimum humidity	5	Maximum humidity	93

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--

08/10/2025 00:51:54 – For information purposes only, can be subject to change