

AUTOTENSIO BPM ARM CUFF ADULT (M/L) BLUE SPENGLER

Ref. SPG440-



CE
0459



08/10/2025 02:51:33 – For information purposes only, can be subject to change

AUTOTENSIO BPM ARM CUFF ADULT (M/L) BLUE SPENGLER



Ref. SPG440-

Ease of use in everyday life for self-measurement family use

Description

AUTOTENSIO® upper arm model. For occasional self-monitoring. - LCD screen, - W.H.O risk indicator- Irregular heart beat detection, - Low batteries indicator, - 2 x 60 memories with date and time, - Purple cuff- Upper arm circumference 22-36 cm size M. Supplied with 4 LR06 AA batteries and a storage bag. Optional power adapter.

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446022501	HS Code	90189010
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	13	Packed unit / Height (cm)	10
Packed unit / Length (cm)	17	Packed unit / Weight	0.594
Bare Unit / Weight (kg)	0.358		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	90

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--