



Ref. 507219



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
5
5 YEARS
GUARANTEE



07/10/2025 23:46:39 – For information purposes only, can be subject to change

Ref. 507219

Spengler takes care of newborns!

Description

The dual Magister Neonatal est dedicated to newborn auscultation.

With its dual head chestpiece equipped of a floating diaphragm on the large side, Magister Neonatal is perfectly adapted to babies for listening high and low frequencies.

Its PVC non-chill rings on both sides of the chestpiece allows a perfect comfort for the newborns.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446042325	HS Code	90181990
Warranty (months)	60		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	16	Packed unit / Height (cm)	4
Packed unit / Length (cm)	32	Bare Unit / Length (cm)	76
Packed unit / Weight	0.25		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	40
Minimum humidity	NA	Maximum humidity	95

OTHER FEATURES

07/10/2025 23:46:39 – For information purposes only, can be subject to change