


DDM ELECTRONIC PNEUMATIC TOURNIQUET EASY PUMP® DUAL 2 PRESSURE CIRCUITS

Ref. G10904



CE
0459



08/10/2025 02:32:02 – For information purposes only, can be subject to change

DDM ELECTRONIC PNEUMATIC TOURNIQUET EASY PUMP® DUAL 2 PRESSURE CIRCUITS

Ref. G10904

The only electronic tourniquet manufactured in France!

Description

The EASY PUMP® DUAL tourniquet is DDM's new generation of electronic tourniquet providing a safe, bloodless operating field for lower limb surgery. A model with 2 independent pneumatic circuits to be used with a single or double cuff or with 2 single cuffs. It allows 2 cuffs to be inflated at different pressures for bilateral surgery or IntraVeinous Regional Anesthesia (IVRA).

Easy to use and intuitive:

- Pre-settings of time, pressure, languages...
- Large screen 7" with IPS touch screen technology (1024 x 600 pixels) with no viewing angle
- Intense brightness (700cd/m²)
- Highly intuitive icons for setup and use.
- Distinct color codes for different circuits

Safe for the patient:

- Integrated LOP function without SpO2 sensor, with corrective ajustement.
- Autotest upon switch on
- Time and pressure preset by default and can be modified during operation.
- IVRA mode and possible guidance function.
- Leakage compensation.
- Visual and audible alarms.

DDM ELECTRONIC PNEUMATIC TOURNIQUET EASY PUMP® DUAL 2 PRESSURE CIRCUITS

Ref. G10904

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700468214960	HS Code	90189084
Warranty (months)	24	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	19
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	25
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	30
Packed unit / Weight	nc	Bare Unit / Weight (kg)	5.5

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	60
Minimum humidity	20	Maximum humidity	80

OTHER FEATURES