

TEMPO MAX HYPOTHERMIC DIGITAL THERMOMETER SPENGLER

Ref. THE002



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2
2 YEARS
GUARANTEE

Description

08/10/2025 02:49:47 – For information purposes only, can be subject to change

TEMPO MAX HYPOTHERMIC DIGITAL THERMOMETER SPENGLER

Ref. [THE002](#)

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	6970392219271	HS Code	90251920
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	2.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	14.3
Packed unit / Weight	nc	Bare Unit / Weight (kg)	0.02

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	95

OTHER FEATURES