

CA-MI NEW ASKIR SUCTION UNIT 30 POLYCARBONATE JAR 1L

Ref. BCASPI12

CA-MI



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2 ANS
2 YEARS
GUARANTEE

Fixed electric surgical Suction Unit, intended to drain body fluids and to perform oral, nasal and tracheal aspiration in adults and children.

Description

08/10/2025 19:51:10 – For information purposes only, can be subject to change

CA-MI NEW ASKIR SUCTION UNIT 30 POLYCARBONATE JAR 1L

CA-MI

Ref. BCASPI12

Ca-mi New Askir 30 Suction Unit, autoclavable polycarbonate jar, 1000 ml

Fixed electric surgical Suction Unit, intended to drain body fluids and to perform oral, nasal and tracheal aspiration in adults and children.

Suitable for tracheostomy patients and is ideal for minor surgical applications as well as postoperative home treatments.

Main applications: Endoscopy, Gastroenterology, General surgery, Home care, Obstetrics, Room, Veterinary Med.

-

The main system is equipped with a pressure gauge (bar & kPa) and a vacuum regulator

-

Optional roller foot-

Supplied accessories: Fluid collection jar with an anti-reflux valve with different options / Silicone hoses (autoclavable) 6 x 10mm and conical connectors ø 8-9-10 mm / Antibacterial & hydrophobic filter (single patient) / CH20 cannula / Power cord with Schuko plug

- Engine: Oil-free piston pump requiring no maintenance

-

Power supply: 230V-50Hz - other configurations available on request

-

Maximum vacuum (adjustable): 0.80 Bar -80 kPa -600 mmHg

-

Max free air flow: 40 l/min

-

Noise level: 60.5 dBA

- Uninterrupted operating cycle

- Weight: 3.6 kg / Dimensions: 35 x 21 x 18 cm

CA-MI NEW ASKIR SUCTION UNIT 30 POLYCARBONATE JAR 1L

CA-MI

Ref. BCASPI12

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	8033622540120	HS Code	90189084
Warranty (months)	24	Country	Italy

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	3.8
Bare Unit / Weight (kg)	3.6		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	45
Minimum humidity	1	Maximum humidity	85

OTHER FEATURES