

AUDIOMETER SIBELMED SIBELSOUND DUO A OPTION AUTO HUGHSON WESTLAKE

Ref. DSI009133

Sibelmed



CE
0318

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
1 AN
YEAR
GUARANTEE



07/10/2025 23:42:24 – For information purposes only, can be subject to change

**AUDIOMETER SIBELMED SIBELSOUND DUO
A OPTION AUTO HUGHSON WESTLAKE****Ref. DSI009133**

A new concept of audiometry for occupational medicine, general care and schools.

Description

Portable and lightweight, the SIBELSOUND DUO audiometer is the ideal device for performing audiometry tests in any environment.

Performs audiometry tests in automatic mode, which requires minimal technician intervention.

Follows the Hughson-Westlake automatic audiometry protocol.

Optional: export test files in standard formats (PDF, XML and CSV).

Features

CE Marking	CE0318	Medical device	Class IIa
Bar code	3661041010754	HS Code	90181990
Warranty (months)	12	Country	Spain

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	32	Packed unit / Height (cm)	14
Packed unit / Length (cm)	44	Packed unit / Weight	2.7
Bare Unit / Weight (kg)	0.6		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES