

Ref. E1.342.012

LUXAMED®



CE



Description

Available in size 2 (E1.342.012), size 3 (E1.343.012) and size 4 (E1.344.012)

08/10/2025 00:36:10 – For information purposes only, can be subject to change

Ref. E1.342.012

LUXAMED®

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381963394	HS Code	90189084
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Size	Size 2		

08/10/2025 00:36:10 – For information purposes only, can be subject to change