


EMERAIR OXYGEN EMERGENCY BAG, RED

Ref. MEEME03_R

ELITE BAGS



07/10/2025 23:10:03 – For information purposes only, can be subject to change

EMERAIR OXYGEN EMERGENCY BAG, RED

Ref. MEEME03_R



For emergency services specialized in first aid

Description

- Interior totally adaptable to the needs of the user thanks to extractable velcro separators
- Elastic bands in all interior flat surfaces
- Reflective strips night and day- 2 external pockets of medium size, 2 external pockets of large capacity, with elastic bands
- 4 extractable compartments with the transparent front part differentiated according to the standardized color code (blue: respiratory system, red: circulatory system, green: various surgical equipment, yellow: pediatric equipment)
- Isothermal amp maker with 27 elastics
- Shoulder strap, handles that surround the bag for better resistance, concealed backpack straps and banana
- Material: 1000D polyamide
- Weight: 3.65 kg- Outdoor measures: 30 x 50 x 23 cm / 32.96 L
- Recommended max. capacity: 25 kg
- Max. capacity of an oxygen cylinder 3L
- Delivered with an ampoulier
- Optional trolley
- Available in waterproof version
- Sold empty

EMERAIR OXYGEN EMERGENCY BAG, RED



Ref. MEEME03_R

Features

CE Marking	NONE	Medical device	No
Bar code	8436564540133	HS Code	4202929890
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	32	Bare Unit / Width (cm)	30
Packed unit / Height (cm)	19	Bare Unit / Height (cm)	23
Packed unit / Length (cm)	50	Bare Unit / Length (cm)	50
Packed unit / Weight	4.1	Bare Unit / Weight (kg)	3.65

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---