


## MEDIC'S ELITE BAGS, BLUE

Ref. MEMEDI02\_B

**ELITE BAGS**



07/10/2025 22:54:53 – For information purposes only, can be subject to change

## MEDIC'S ELITE BAGS, BLUE

Ref. MEMEDI02\_B



### Description

Sanitary sports case with practical and comfortable design. Its internal organization totally adaptable thanks to the extractable velcro dividers, allows the rapid location of the equipment in case of medical or sports emergency.

- Quilted handle with velcro closure- Card holder
- Reinforced structure with ribbon
- Isothermal amplestick with 41 elastic strips- Document holder
- Large capacity external pocket with elastic bands and spacious compartments
- Adjustable honeycomb structure interior
- Flat pocket for the document holder
- 2 flat net pockets for instruments
- External pocket for personal oblets
- Measures: 50 x 25 x 31 cm
- Capacity: 33.6L
- Weight: 2.23 kg
- Material: 1200D polyester. Washable
- Recommended max. capacity: 12 kg
- Optional trolley
- Sold empty
- Available in blue (MEMEDI02\_B) and black (MEMEDI02\_N)

## MEDIC'S ELITE BAGS, BLUE

Ref. MEMEDI02\_B



### Features

CE Marking	NONE	Medical device	No
Bar code	8436564540423	HS Code	42029270
Warranty (months)	12	Country	China

#### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	25
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	31
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	50
Packed unit / Weight	3.55	Bare Unit / Weight (kg)	2.23

#### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

#### OTHER FEATURES

Colors	Blue	Packed unit / Quantity	1
--------	------	------------------------	---