

PLASTIC TROLLEY, 2 TRAYS 61,5 X 45 CM

Ref. OM614501



PLASTIC TROLLEY, 2 TRAYS 61,5 X 45 CM

Ref. [OM614501](#)

Lightweight, handy and very robust pedestal table

Description

The quality of the casters and ergonomic handles ensure maximum comfort when travelling.

- External dimensions: 80 x 48 x H 95 cm, including handles

- Dimensions of the trays: 61.5 x 45 cm

-Maximum load: 180 kg per trolley and 60 kg per level

-Color: Light gray

-Materials: Aluminum upright, PP top

- 4 PVC swivel casters, Ø 100 mm including 2 brake

-Weight: 8.5 kg

Features

CE Marking	NONE	Medical device	No
Bar code	3760206421227	HS Code	39269097
Warranty (months)	12		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	71	Bare Unit / Width (cm)	48
Packed unit / Height (cm)	19	Bare Unit / Height (cm)	95
Packed unit / Length (cm)	52	Bare Unit / Length (cm)	80
Packed unit / Weight	9.806	Bare Unit / Weight (kg)	8.5

SPECIFIC STORAGE CONDITIONS

Minimum temperature	5	Max temperature	35
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---