

ALUMINIUM FOLDABLE STRETCHER IN 2 (LENGTH AND LARGE) + POUCH, ORANGE

Ref. SBRAN08_O



ALUMINIUM FOLDABLE STRETCHER IN 2 (LENGTH AND LARGE) + POUCH, ORANGE

Ref. SBRAN08_O



Folding stretcher with balance wheel system

Description

Folding stretcher with balancing system

- High resistance aluminium structure and waterproof fabric, fire resistant and anti-cracking
- This balancing stretcher folds in length and width
- Lightweight, easy to carry
- Dimensions (L x W x H): 208 X 55 x 13 cm, unfolded;
104 x 17 x 9 cm, folded
- Weight: 8 kg
- Maximum load carried: 159 kg
- Single colour: orange
- Delivered with a stretcher bag

ALUMINIUM FOLDABLE STRETCHER IN 2 (LENGHT AND LARG) + POUCH, ORANGE

Ref. SBRAN08_O



Features

CE Marking	CE	Medical device	Class I
Bar code	3661040094632	HS Code	94029000
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	20	Bare Unit / Width (cm)	54.5
Packed unit / Height (cm)	10	Bare Unit / Height (cm)	13.5
Packed unit / Length (cm)	10	Bare Unit / Length (cm)	208
Packed unit / Weight	5.6	Bare Unit / Weight (kg)	4.7

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---