

MALE CONNECTOR (FOR VACUUM PORT) FOR ASKIR 20 PACK OF 5'S

CA-MI

Ref. BCASPI164



CE



📅 12/10/2025 12:19:17 – For information purposes only, can be subject to change

MALE CONNECTOR (FOR VACUUM PORT) FOR ASKIR 20 PACK OF 5'S



Ref. [BCASPI164](#)

Features

CE Marking	CE	Medical device	No
Bar code	3661040160375	HS Code	90330090
Warranty (months)	12	Country	Italy

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.01
Bare Unit / Weight (kg)	0.01		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---