

Ref. BT103S

CE



12/10/2025 12:19:51 – For information purposes only, can be subject to change

Ref. BT103S

Description

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468215684	HS Code	9018901000
Country	France		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES