

CUSHING RETRACTOR 25CM 30MM

Ref. IE08830



Features

CE Marking	CE	Medical device	Class I
Bar code	3661040016443	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	2.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	24.6
Packed unit / Weight	0.059	Bare Unit / Weight (kg)	0.059

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---