

ELITE BAGS - MEDICAL LEATHER BRIEFCASE CLASSY'S

Ref. MECLASS01_B

 **ELITE BAGS**



 12/10/2025 13:07:13 – For information purposes only, can be subject to change

ELITE BAGS - MEDICAL LEATHER BRIEFCASE CLASSY'S



Ref. [MECLASS01_B](#)

Medical leather bag with traditional design.

Description

With compact size to carry essential diagnostic equipment. A professional touch of elegance and distinction.

- Internal pocket for personal items
- With key clasp
- With handle and shoulder strap
- Material: leather
- Outdoor measures: 36 x 21 x 20 cm / 15.12 L
- Weight: 1.15 kg
- Recommended max. capacity: 3 kg
- Sold empty
- Available in brown (MECLASS01) and blue (MECLASS01)

Features

CE Marking	NONE	Medical device	No
Bar code	8436564542571	HS Code	42029900
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	21
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	20
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	36
Packed unit / Weight	1.75	Bare Unit / Weight (kg)	1.15

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Blue	Packed unit / Quantity	1
--------	------	------------------------	---