

## MEDIPOLOGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE

Ref. WSAM\_45L

medi Plogs



📅 12/10/2025 11:28:08 – For information purposes only, can be subject to change

**MEDIPOLOGS OPERATING ROOM SHOES,  
NON-PERFORATED, 35, NAVY BLUE**Ref. WSAM\_45L

=====

## Description

Clogs suitable for modern hospitals, operating theatres and intensive care units.

Products designed and manufactured in collaboration with researchers and doctors.

- Antistatic clogs with non-slip and non-chafing soles
- The back of the heel is sloped and rounded to make walking easier
- Ventilation system allows air to circulate
- Toe prints to enhance hold and spurs to aid the circulation
- 100% CFC-free polyurethane
- Washable by hand or machine at 60°C
- Sizes: 35 - 47

## MEDILOGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE



Ref. WSAM\_45L

### Features

|                             |               |                           |          |
|-----------------------------|---------------|---------------------------|----------|
| CE Marking                  | NONE          | Medical device            | No       |
| Bar code                    | 3661040100319 | HS Code                   | 64029993 |
| Country                     | Germany       |                           |          |
| DIMENSIONS AND WEIGHT       |               |                           |          |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 1.16     |
| Bare Unit / Weight (kg)     | 1             |                           |          |
| SPECIFIC STORAGE CONDITIONS |               |                           |          |
| Minimum temperature         | NA            | Max temperature           | NA       |
| Minimum humidity            | NA            | Maximum humidity          | NA       |
| OTHER FEATURES              |               |                           |          |
| Colors                      | fuchsia       | Shoe Size                 | 45       |