

Ref. BM101S



CE



Description

14/10/2025 23:49:59 – For information purposes only, can be subject to change

Ref. **BM101S**

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468215721	HS Code	9018901000
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100
OTHER FEATURES			