

Ref. **BM105S**

Description

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468215769	HS Code	9018901000
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100
OTHER FEATURES			