


Ref. TI50003S

Spengler



CE
0459



14/10/2025 23:45:12 – For information purposes only, can be subject to change

Ref. TI50003S



Description

Features

CE Marking	CE0459	Medical device	No
Bar code	3700468215516	HS Code	90189010
Warranty (months)	36	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.485
Bare Unit / Weight (kg)	0.325		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES

14/10/2025 23:45:12 – For information purposes only, can be subject to change