

Ref. TI50004S

 Spengler



CE
0459



 14/10/2025 23:44:09 – For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr

Ref. TI50004S



Description

Features

CE Marking	CE0459	Medical device	No
Bar code	3700468215523	HS Code	90189010
Warranty (months)	36	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES

14/10/2025 23:44:09 – For information purposes only, can be subject to change