

Ref. TI50005S

Spengler



CE
0459

GARANTIE
3 ANS
3 YEARS
GUARANTEE



14/10/2025 23:08:11 – For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

Ref. TI50005S



Description

Features

CE Marking	CE0459	Medical device	No
Bar code	3700468215530	HS Code	90189010
Warranty (months)	36	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES

14/10/2025 23:08:11 – For information purposes only, can be subject to change