



Ref. TS53000S



CE  
0459



14/10/2025 23:26:59 – For information purposes only, can be subject to change

Ref. TS53000S



## Description

## Features

|            |               |                |          |
|------------|---------------|----------------|----------|
| CE Marking | CE0459        | Medical device | No       |
| Bar code   | 3700468215486 | HS Code        | 90189010 |
| Country    | France        |                |          |

### DIMENSIONS AND WEIGHT

|                           |    |                           |    |
|---------------------------|----|---------------------------|----|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | nc |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |     |
|---------------------|----|------------------|-----|
| Minimum temperature | 10 | Max temperature  | 70  |
| Minimum humidity    | 10 | Maximum humidity | 100 |

### OTHER FEATURES