

Ref. 335111



Description

Cuffs with dual tube bladder & coupling (60 cm tube with Classic® bulb and valve, and 10 cm with male screw-on metal fitting)

Available in :

Adult (M), nylon, velcro, black (335111)

Adult (L), nylon, black (335101)

18/10/2025 08:25:21 – For information purposes only, can be subject to change

Ref. 335111



Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028435	HS Code	9018901000
Warranty (months)	6	Country	Indonesia

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.245

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---