




LED FIBER OPTIC OTOSCOPE - SMARTLED 5500

Ref. 570517C



20/10/2025 20:24:30 – For information purposes only, can be subject to change

LED FIBER OPTIC OTOSCOPE - SMARTLED 5500

Ref. [570517C](#)



The SMARTLED 5500 otoscope has been developed with a lens to allow exceptional diagnosis.

Description

LED lighting with Optical fibreThe Spengler otoscope SMARTLED 5500® is using LED lighting with optical fibre allowing a direct, homogeneous illumination and without reflections of the auricular canal and the eardrum. Luminance 85K LUX. Bimaterial head in ABS and Stainless steel, screwed head for a safe fixation. Ergonomic ABS grip. End clamp with built-in switch button. Optical fibre for optimal transmission of the light beam. x 3 magnification. Supplied in a case with two LR6 AA batteries, 8 speculums (two of 2.5 mm and four of 4 mm).

Available in 3 colors: black (570517), coral (570517C), blueberry (570517B) and also in a paediatric version (light green 570517P)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446041700	HS Code	90189042
Warranty (months)	36	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	12	Packed unit / Height (cm)	4
Packed unit / Length (cm)	20	Packed unit / Weight	0.385
Bare Unit / Weight (kg)	0.075		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

Colors	Coral	Packed unit / Quantity	1
--------	-------	------------------------	---