

## LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)

Ref. 513210



CE  
0459



**Light and robust, the blood pressure monitor that combines style and performance.**

### Description

MANOMETRE EN CONTROLE QUALITE : VERIFIEZ QUE VOTRE LOT EST ACCEPTE, AVERTIR POUR CONTROLE PAR ECHANTILLONNAGE

21/10/2025 22:37:10 – For information purposes only, can be subject to change

## LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)

Ref. 513210



### Features

|                             |               |                           |                               |
|-----------------------------|---------------|---------------------------|-------------------------------|
| CE Marking                  | CE0459        | Medical device            | Class Im (measuring function) |
| Bar code                    | 3700446009953 | HS Code                   | 90189010                      |
| Warranty (months)           | 24            | Country                   | France                        |
| DIMENSIONS AND WEIGHT       |               |                           |                               |
| Packed unit / Width (cm)    | 13            | Packed unit / Height (cm) | 7                             |
| Packed unit / Length (cm)   | 19            | Packed unit / Weight      | 0.5                           |
| Bare Unit / Weight (kg)     | 0.4           |                           |                               |
| SPECIFIC STORAGE CONDITIONS |               |                           |                               |
| Minimum temperature         | -20           | Max temperature           | 70                            |
| Minimum humidity            | 85            | Maximum humidity          | 85                            |
| OTHER FEATURES              |               |                           |                               |
| Packed unit / Quantity      | 1             |                           |                               |