

LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)

Ref. 513229



CE
0459



Light and robust, the blood pressure monitor that combines style and performance.

Description

MANOMETRE EN CONTROLE QUALITE : VERIFIEZ QUE VOTRE LOT EST ACCEPTE, AVERTIR POUR CONTROLE PAR ECHANTILLONNAGE

22/10/2025 00:21:15 – For information purposes only, can be subject to change

LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)

Ref. [513229](#)



Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446007195	HS Code	90189010
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	7
Packed unit / Length (cm)	19	Packed unit / Weight	0.526
Bare Unit / Weight (kg)	0.5		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		