

LIAN METAL / GREY COTTON CUFF / ADULT (M)

Ref. 521318



CE
0459



21/10/2025 21:56:15 – For information purposes only, can be subject to change

**Ref. 521318**

Ambidextrous model offering great comfort of use thanks to its large dial.

Description

LIAN METAL Velcro armband (XXS to XL) Spengler
Ambidextrous model offering great comfort of use thanks to its large dial.

>>> WARNING:
Possible break in some sizes, do not hesitate to consult the sales department for equivalencies.

HIGH PRECISION OF MEASUREMENTS• High precision "watchmaking" mechanism

COMFORT• Large dial for simplified reading•
Ergonomic spoon and grooved pear for a perfect grip•
Ambidextrous by simple rotation•
Large volume pear for fast inflation
• Velcro Spengler closure, the perfect mix between support and practicality

**SHOCK
RESISTANT**• Non-deformable metal housing•
Shockproof foam pressure gauge protector

• Available in Clinic® version for prevention of cross-contamination

LIAN METAL / GREY COTTON CUFF / ADULT (M)

Ref. 521318



Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446008840	HS Code	90189010
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.75
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		