

SET OF 20 SUNTECH MONOBLOC
ARMBANDS DOUBLE TUBING SINGLE
PATIENT SINGLE ADULT LONG (23-33CM) -
SCREW CONNECTOR

SUNTECH

Ref. DSU231086



CE



**SET OF 20 SUNTECH MONOBLOC
ARMBANDS DOUBLE TUBING SINGLE
PATIENT SINGLE ADULT LONG (23-33CM) -
SCREW CONNECTOR**

SUNTECH

Ref. DSU231086

Pack of 20 Suntech single-piece double-tube adult cuffs – precision, safety, and comfort for every unique patient.

Description

Set of 20 Suntech "SOFT" single-patient armbands, royal blue.
Double tubing, adult size 23-33 cm long, screw connection.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446014537	HS Code	9018905040
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.91

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---