

IDEAL+ SINGLE HEAD ADULT, BLUE

Ref. DH7081_N



IDEAL+ SINGLE HEAD ADULT, BLUE

Ref. DH7081_N



Single-head adult stethoscope, designed for reliable and accurate auscultation.

Description

- Stethoscope with single adult metal horn (Ø: 45 mm)•
- Anti-cold ring and horn of the color of the tube•
- Aluminum lyre incorporated into the tubing
- Also available in double pavilion for adults and a pediatric model

Available in Black (DH7081_N) or Red (DH7081_R)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661041000120	HS Code	9018908400
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.15
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
OTHER FEATURES			
Colors	Black	Packed unit / Quantity	1