

MONOBLOC CONICAL ADULT ARM CUFF WITH ONE BLADDER FOR DDM TOURNIQUETS



Ref. GBCM105

The monobloc conical arm cuff with one bladder!

Description

This **Adult L** monobloc arm cuff made of **PU** is decontaminable. With a **conical design** and only **one bladder** it is made c a monobloc unit with a ribbon to secure the attachment.

Width of the bladder : 10cm

Length of the bladder : 47cm

Total length : 70cm

The arm cuff can not be sterilized. It can be dry-cleaned or immersed with paying attention to maintain he connector of the tube out of water.

It is also available as lower limb cuffs

**MONOBLOC CONICAL ADULT ARM CUFF
WITH ONE BLADDER FOR DDM
TOURNIQUETS****Ref. GBCM105****Features**

CE Marking	CE	Medical device	Class I
Bar code	3700468206699	HS Code	90189084
Country	France		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES