

CONICAL LOWER LIMB CUFF WITH 1 BLADDER FOR DDM TOURNIQUET

Ref. GCCS105



Those cuffs are autoclavable!

Description

Cuffs and lower limb cuff for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of **1 removable silicon bladder** for easy autoclaving.

Those ones are designed in a **conical form**

Several sizes are available

Dimensions of the bladder in cm

Lower limb cuff adult (GCCS104)	69 x 9 cm
Lower limb cuff adult L (GCCS105)	81 x 9 cm
Lower limb cuff adult XL (GCCS106)	85,5 x 12 cm

Dimensions de la poche en cm

Cuissard adulte (GCCS104)	69 x 9 cm
Cuissard adulte L (GCCS105)	81 x 9 cm
Cuissard adulte XL (GCCS106)	85,5 x 12 cm

**CONICAL LOWER LIMB CUFF WITH 1
BLADDER FOR DDM TOURNIQUET**

Ref. GCCS105

**Features**

CE Marking	CE	Medical device	Class I
Bar code	3700468203568	HS Code	90189084
Warranty (months)	24	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES