

CHILD TIPS FOR PEAK FLOWMETERS Ø 22 MM - BOX OF 100 PCS

Ref. L35678



CE



24/10/2025 21:49:27 – For information purposes only, can be subject to change

CHILD TIPS FOR PEAK FLOWMETERS Ø 22 MM - BOX OF 100 PCS

Ref. L35678



Box of 100 child mouthpieces (Ø22 mm) for Peak Flow meter, designed for hygienic and safe use.

Description

Pediatric tips for SPENGLER peak flow meter.
Dimensions of the tips 22mm. Length : 67 mm
Sold in boxes of 100 pieces.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446001698	HS Code	48229000
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.584
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		