

SPO2 SENSOR PEDIATRIC FINGER CLIP FOR MASTER PALM 2

Ref. 222251



CE
0537

DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
1 AN
1 YEAR
GUARANTEE



27/10/2025 21:34:30 – For information purposes only, can be subject to change

SPO2 SENSOR PEDIATRIC FINGER CLIP FOR MASTER PALM 2



Ref. 222251

Used for oxygen saturation measurement in children.

Description

Rigid pediatric SpO2 sensor for MASTER PALM 2 oximeter

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446030216	HS Code	90181990
Warranty (months)	12	Country	Uruguay
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.093
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	1	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 27/10/2025 21:34:30 – For information purposes only, can be subject to change