



## SPO2 SENSOR PEDIATRIC SILICON FOR MASTER PALM 2

Ref. 222255



CE  
0537



📅 27/10/2025 21:32:55 – For information purposes only, can be subject to change

## SPO2 SENSOR PEDIATRIC SILICON FOR MASTER PALM 2



Ref. 222255

Used for oxygen saturation measurement in children.

### Description

Pediatric silicon spo2 sensor for oximeter MASTER PALM 2

### Features

|                   |               |                |           |
|-------------------|---------------|----------------|-----------|
| CE Marking        | CE0537        | Medical device | Class IIb |
| Bar code          | 3700446030254 | HS Code        | 90181990  |
| Warranty (months) | 12            | Country        | Uruguay   |

#### DIMENSIONS AND WEIGHT

|                           |    |                           |       |
|---------------------------|----|---------------------------|-------|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc    |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | 0.084 |

#### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | 1  | Max temperature  | 50 |
| Minimum humidity    | NA | Maximum humidity | NA |

#### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|