

PAIR OF CORAL EARTIPS FOR PULSE II

Ref. 506903



Designed for otoscopic use.

Description

Pair of ear tips for SPENGLER PULSE II. Standard size, coral color.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030018	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.004
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		