

BURGUNDY CHROME PLATED DUAL TUBE BINAURALS



Ref. 507923

burgundy chrome plated dual tube binaurals

Description

Double-tube chrome lyre in burgundy, designed by Spengler for optimal acoustic transmission and superior listening comfort. Perfect for demanding healthcare professionals.

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446009601	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.08
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		