

FLOATING DIAPHRAGM FOR MAGISTER II SPENGLER

Ref. 507975



Used in manual blood pressure measurement.

Description

Spare part for stethoscopes MAGISTER® II from SPENGLER.
Multi-frequency fiberglass membrane for adult pavilion.

📅 27/10/2025 21:28:20 – For information purposes only, can be subject to change

FLOATING DIAPHRAGM FOR MAGISTER II SPENGLER



Ref. 507975

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030346	HS Code	90181990
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.001
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		