

EXTRA LARGE ADULT CUFF 40-48CM FOR SPOT-CHECK 2VS

Ref. 546021



Used for non-invasive blood pressure measurement.

Description

Extra large armbands (40-48cm) for SPENGLER 2VS spot-check station.

27/10/2025 21:41:31 – For information purposes only, can be subject to change

EXTRA LARGE ADULT CUFF 40-48CM FOR SPOT-CHECK 2VS

Ref. 546021



Features

CE Marking	CE	Medical device	Class I
Bar code	3700446029685	HS Code	9018901000
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.102

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---