

SET OF 4 PAPER ROLLS FOR SPENGLER 4VS MONITOR PRINTER

Ref. 551643



Used for printing monitoring data.

Description

Set of 4 rolls of SPENGLER 4VS Vital Signs Monitor Printer Paper

27/10/2025 21:29:40 – For information purposes only, can be subject to change

SET OF 4 PAPER ROLLS FOR SPENGLER 4VS MONITOR PRINTER



Ref. 551643

Features

CE Marking	NONE	Medical device	No
Bar code	3700446033125	HS Code	85049099
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.228
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		