

REPLACEMENT HANDLE FOR SMARTLED 5500 AND VISIONLED 4000

Ref. 570519



Intended for medical lighting use.

Description

Replacement handle for otoscope SPENGLER SMARTLED 5500

 27/10/2025 21:37:08 – For information purposes only, can be subject to change

REPLACEMENT HANDLE FOR SMARTLED 5500 AND VISIONLED 4000



Ref. 570519

Features

CE Marking	NONE	Medical device	No
Bar code	3700446026370	HS Code	90189042
Country	Pakistan		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.038
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		