

SET OF 5 SUNTECH MONOBLOC ARMBANDS REUSABLE CHILD (12-19CM) - SCREW CONNECTOR

SUNTECH

Ref. [DSU231001](#)



CE



Enables oxygen saturation monitoring in children.

Description

Reusable nylon cuff sold by 5. Single tube, screw-on coupling.

27/10/2025 21:34:00 – For information purposes only, can be subject to change

SET OF 5 SUNTECH MONOBLOC ARMBANDS REUSABLE CHILD (12-19CM) - SCREW CONNECTOR

SUNTECH

Ref. DSU231001

Features

CE Marking	CE	Medical device	Class I
Bar code	0120840935100383	HS Code	9018905040
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.216

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---