

LOT 5 SUNTECH LAR.ADUL NYLON ARMBANDS (BURGUNDY) 31-40CM WITHOUT RACC

Ref. [DSU231027](#)

SUNTECH



CE



Intended for blood pressure measurement.

Description

Reusable nylon cuff sold by 5. Single tube, no coupling.

27/10/2025 21:26:29 – For information purposes only, can be subject to change

**LOT 5 SUNTECH LAR.ADUL NYLON
ARMBANDS (BURGUNDY) 31-40CM
WITHOUT RACC**

SUNTECH

Ref. DSU231027

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446013769	HS Code	9018905040
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	23	Packed unit / Height (cm)	10
Packed unit / Length (cm)	25	Packed unit / Weight	0.616
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		