

SET OF 20 SUNTECH MONOBLOC ARMBANDS SINGLE PATIENT LONG CHILD (12-19CM) - BAYONET CONNECTOR

SUNTECH

Ref. [DSU231062](#)

Designed for occasional pediatric blood pressure measurement.

Description

Set of 20 Suntech "SOFT" single-patient armbands, green.
Child size long 12-19cm, rectus bayonnette connection.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446014100	HS Code	9018905040
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.857

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---