

LARGE ADULT 31-40CM X20

Ref. DSU231087

SUNTECH



CE



Intended for single-use blood pressure measurement in adults.

Description

Set of 20 Suntech "SOFT" single-patient armbands, burgundy.
Double tubing, Wide adult size 31-40 cm, screw connection.

27/10/2025 21:35:29 – For information purposes only, can be subject to change

LARGE ADULT 31-40CM X20

Ref. [DSU231087](#)

SUNTECH

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446014544	HS Code	9018905040
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	2.145
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		