

SET OF 5 SUNTECH MONOBLOC REUSABLE ARMBANDS EXTRA LARGE (38-50CM) - BAYONET METAL CONNECTOR

SUNTECH

Ref. [DSU231116](#)



CE



Intended for non-invasive blood pressure measurement in adults.

Description

Reusable nylon cuff sold in boxes of 5 units.
Single tubing, bayonette connection 38-50cm.

27/10/2025 21:22:53 – For information purposes only, can be subject to change

SET OF 5 SUNTECH MONOBLOC REUSABLE ARMBANDS EXTRA LARGE (38-50CM) - BAYONET METAL CONNECTOR

SUNTECH

Ref. [DSU231116](#)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661041008843	HS Code	9018905040
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	25	Packed unit / Height (cm)	10
Packed unit / Length (cm)	23	Packed unit / Weight	0.8

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---