

10'/3M PATIENT BP HOSE BAYONET TO BAYONET

SUNTECH

Ref. [DSU551456](#)

Used for connecting blood pressure cuffs.

Description

Pediatric silicone SpO₂ sensor compatible with the Spengler Master Palm 2 monitor. Provides comfortable and reliable oxygen saturation measurement in children.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446032562	HS Code	3917.23.0000
Warranty (months)	12		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	2	Packed unit / Height (cm)	18
Packed unit / Length (cm)	19	Packed unit / Weight	0.1

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---