



SPO2 SENSOR NEONATE Y FOR MASTER PALM 2

Ref. 222252



CE
0537



29/10/2025 02:13:39 – For information purposes only, can be subject to change

SPO2 SENSOR NEONATE Y FOR MASTER PALM 2

Ref. 222252



A Y-shaped neonatal sensor designed for gentle and accurate SpO₂ monitoring with the Master Palm 2.

Description

SpO₂ neonate y sensor for oximeter MASTER PALM 2

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446030223	HS Code	90181990
Warranty (months)	12	Country	Uruguay

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.086

SPECIFIC STORAGE CONDITIONS

Minimum temperature	1	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---