


SPO2 SENSOR ADULT SILICON FOR MASTER PALM 2

Ref. 222254



CE
0537



29/10/2025 02:17:13 – For information purposes only, can be subject to change

SPO2 SENSOR ADULT SILICON FOR MASTER PALM 2

Ref. 222254



A soft and accurate sensor for reliable daily SpO₂ monitoring with optimal comfort.

Description

Adult silicone spo2 sensor for oximeter MASTER PALM 2

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446030247	HS Code	90181990
Warranty (months)	12	Country	Uruguay

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.099

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---